

OPTIMAL SAMPLE CONDITIONS FOR ApoE QUANTIFICATION BY APO E (Spinreact)

| | |
|-----------------------------|--------------------------------------|
| Matrix: | EDTA plasma EDTA heparin Serum |
| Recommended minimum volume: | 200 µL |
| Storage conditions: | Frozen (3 months at -20°C) |

If samples do not match these conditions, please contact services@biocross.es.

APO E is an in vitro diagnostic test for quantification of the Apolipoprotein E (ApoE). APO E is CE-marked and is commercially distributed by Spinreact. Biocross is providing this analysis only for Research purposes.